



Faulding® ReduZnore Mouth Piece 1 Unit

Faulding® ReduZnore Mouth Piece A mouldable mouthpiece designed to help relieve snoring also:

- Helps relieve snoring.
- Improves air flow.
- Widens nasal & throat passages.
- Compatible iPhone app available for purchase.
- Bonus nasal inserts included.

About this product

Faulding® ReduZnore Snoring Mouthpiece is an innovative approach to snoring control. It contains a customizable mouthpiece, complemented by bonus nasal dilators. The mouthpiece is moulded to the users teeth, and once moulded properly, creates a mandibular advancement of between 2-7mm from its normal position.

Creating a mandibular advancement, or positioning the lower jaw slightly forward, helps to:

- Widen the throat, allowing for an increase in air flow
- Put more tension on the soft, loose tissue located at the very back of the mouth & throat which, when it vibrates, contributes to snoring
- Make breathing more calm and relaxed
- Decrease snoring

For best results, its recommended to incorporate the Faulding® MegaVent Nasal Dilators with the Faulding® ReduZnore Snoring Mouthpiece when looking to alleviating snoring. Faulding® MegaVent is a snoring and breathing aid which encourages users to breathe through their nose, allowing for increased nasal air flow. By combining both the ReduZnore Mouthpiece and the MegaVent Nasal Dilators, you will be able to breathe with a calm, continuous airflow.

Care Instructions: Rinse thoroughly in hot tap water before and after use. Soak device once a week in mouthwash for continued freshness. Replace annually.

Compatible iPhone App Available The Faulding® Snoring Monitor iPhone App is a professional snoring aid that includes a sleep lab test providing audio and graphed results of breathing, snoring and chest movements. It can help detect signs of sleep apnoea. It's available for purchase from the iTunes store, it can be found in the 'medical' category or by clicking here: <https://itunes.apple.com/au/app/faulding-snoring-monitor/id641485485?mt=8>



Ingredients

EVA Thermoplastic.

Free from

Latex, silicon.

Contains

EVA Thermoplast (both mouthpiece and dilators).

Directions for use

Follow usage instructions provided with unit.

Warnings

Always read the label. Follow fitting and usage instructions carefully. Not recommended for people with dentures, joint/jaw problems, a large overbite, general teeth problems or people 18 years of age or under. Full instructions available inside the Faulding® ReduZnore Snoring Mouthpiece pack. Use only as directed.

Product Storage

Store below 30°C.

Barcode

9316100686718

Product Dimensions

135mm H x 80mm W x 60mm D